

## PRO RIDER C2 DELUXE

Mobility Scooter available in Blue



- ◆ **AMAZING 72AH TOTAL BATTERY POWER**
- ◆ **FREE 5 STAR HOME WARRANTY COVER**
- ◆ **22 STONE COMFORTABLE USER WEIGHT**
- ◆ **22 MILES RANGE ON ONE CHARGE**



### 4 Mph Class 2 Mobility Scooter

With the Deluxe Scooter you have the absolute freedom of using a mobility scooter which is strong, stable and purpose built for long journeys.

This top of the range mobility scooter has a built in Head Lights and side reflectors, a Deluxe Swivel seat with side adjusters, Extra back rest support for the seat and neck, Height adjustable seat, angle fold down adjustable seat, foot rests, steering wheel angle adjustment, front bumper and detachable basket.

The Pro Rider Deluxe Mobility scooter incorporates the latest in Electromagnetic braking systems, giving the ultimate in user safety.

With it's powerful 560 watt motor and heavy duty 2 x 36ah batteries the Pro Rider Deluxe Mobility scooter can handle just about any incline you will come across.

It's superbly designed frame and seat system as well as the heavy duty tyres will allow the user to navigate through bumps and gradients on the path and roads with ease.

Simple to use fingertip controls make the whole experience of using a mobility scooter both enjoyable and safe for the user. The unique head light system allows for use during darkness.

The new Deluxe Mobility Scooter is a revolutionary design from the worlds leading scooter manufacturer. We believe it is the most user-friendly Top of the range mobility scooter on the market!