

PRO RIDER C3 DELUXE

**Professional Road Mobility Scooter
available in Red**



- ◆ **AMAZING 72AH TOTAL BATTERY POWER**
- ◆ **FREE 5 STAR HOME WARRANTY COVER**
- ◆ **24 STONE COMFORTABLE USER WEIGHT**
- ◆ **24 MILES RANGE ON ONE CHARGE**



6 Mph Class 3 Mobility Scooter

We are proud to present the ultimate in Deluxe Mobility scooter use.

With the C3 Deluxe Scooter you have the absolute freedom of using a mobility scooter which is strong, stable and purpose built for long journeys. This top of the range mobility scooter has built in Head Lights, Rear Lights, Indicators, Hazard Lights and side reflectors, a Deluxe Swivel seat with side adjusters and length adjusters, Extra back rest support for the seat and neck, Height adjustable seat, angle fold down adjustable seat, foot rests, steering wheel angle adjustment, front bumper and detachable basket.

The scooter has a top speed of 6 mph and a Hi Low adjustable switch to enable both pavement and road use.

The Pro Rider C3 Deluxe Mobility scooter incorporates the latest in Electromagnetic braking systems, giving the ultimate in user safety.

With its powerful 650 watt motor and heavy duty 2 x 36ah batteries the Pro Rider C3 Deluxe Mobility scooter can handle just about any incline you will come across. It's superbly designed frame and seat system as well as the heavy duty tyres will allow the user to navigate through bumps and gradients on the path and roads with ease.

Simple to use fingertip controls make the whole experience of using a mobility scooter both enjoyable and safe for the user. The full lights package allows for use during darkness.

The C3 Deluxe Mobility Scooter will provide many years of mobility for you.