

PRO RIDER ELITE

FEATURES

- ◆ 20 Stones max user weight
- ◆ Easy to Dismantle and re-assemble in just seconds
- ◆ Comfortable DELUX extra padded height adjustable seat
- ◆ Width adjustable armrests
- ◆ Fold down arms and backrest
- ◆ Detachable Basket
- ◆ Horn
- ◆ Battery charge indicator
- ◆ Anti tippers
- ◆ Puncture proof tyres
- ◆ Direct access to battery and free wheel mechanism
- ◆ Five Star warranty and service guarantee



SPECIFICATIONS

- ◆ Overall Length 40
- ◆ Overall Width 20" (Width can be reduced by adjusting the arm rests)
- ◆ Overall Height 33"
- ◆ Ground clearance 3"
- ◆ Seat Width 18"
- ◆ Battery Size 2 x 12 AH Sealed
- ◆ Battery pack weight 15 lbs
- ◆ Range Up to 13 miles *
- ◆ Maximum Speed 4 mph
- ◆ Maximum User Weight 20 stone
- ◆ FREE High Powered Charger