

PRO RIDER ROAD KING

FEATURES

- ◆ 160 kgs max user weight
- ◆ Easy to Dismantle and re-assemble in just seconds
- ◆ Comfortable DELUX extra padded height adjustable seat
- ◆ Width adjustable armrests
- ◆ Fold down arms and backrest with swivel seat
- ◆ Detachable Basket
- ◆ Horn / Head lights / Indicators
- ◆ Battery charge indicator
- ◆ Anti tippers
- ◆ Strong tyres to help prevent punctures
- ◆ Direct access to battery and free wheel mechanism
- ◆ Five Star warranty and service guarantee
- ◆ The front handles can be adjusted to fix in several positions
- ◆ Deluxe seat can be adjusted to move forwards, backwards, up, down and forward tilt
- ◆ This scooter has both front and rear lights



SPECIFICATIONS

- ◆ Overall Length 55" (Length can be shortened by removing the basket and anti tip wheels)
- ◆ Overall Width 25.3" (Width can be reduced by adjusting the arm rests)
- ◆ Overall Height 51"
- ◆ Battery Size 2 x 36 AH Sealed
- ◆ Total weight (inc battery, seat and basket) 96 kgs
- ◆ Range Up to 30 miles *
- ◆ Maximum Speed 12 kph / 8 mph
- ◆ Maximum User Weight 160 kgs (26 stone)