

PRO RIDER ROAD KING

**Professional Road Mobility Scooter
available in both Red and Black**



- ◆ **AMAZING 72AH TOTAL BATTERY POWER**
- ◆ **FREE 5 STAR HOME WARRANTY COVER**
- ◆ **26 STONE COMFORTABLE USER WEIGHT**
- ◆ **30 MILES RANGE ON ONE CHARGE**
- ◆ **FULLY SHROUDED P BAR TILLER**



8 Mph Road Class 3 Scooter

With the Road King Class 3 mobility scooter you have the absolute freedom of using a mobility scooter which is strong stable and purpose built for long journeys. This top of the range road mobility scooter has built in head lights and side reflectors, left and right indicators a deluxe swivel seat with side adjusters, extra back rest support for the seat and neck, height adjustable seat, foot rest, steering wheel angle adjustment, front bumper and attachable basket.

The Pro Rider Road King incorporates the latest in Electromagnetic braking systems, giving the ultimate in user safety.

With it's huge powerful 1100 watt motor and heavy duty 2 x 36ah batteries the Pro Rider Road King can handle the most demanding of inclines. It's superbly designed suspension system and heavy duty tyres will allow the user to navigate through the bumps and gradients on the paths and roads. The Road King has the latest Hi-tech super suspension system allowing the user to enjoy the perfect soft ride - perfect for those long journeys.

Simple to use finger-tip controls make the whole experience of using a mobility scooter both enjoyable and safe for the user. The unique headlight system allows for use during darkness.